

From Educational MS to GC-MS Pairing the IQ-2000™ and GC

Fast Fact

“Add-on” to Education by Turning the IQ-2000™ into an Educational GC-MS

Gas Chromatography Mass Spectrometry (GC-MS) is a widely used analytical technique which harnesses the power of two analytical methods. Gas Chromatography separates vaporized compounds without decomposition. Mass spectrometry ionizes compounds and separates them based on a mass-to-charge ratio producing a mass spectrum which can be used for chemical identification and concentration calculations. The combination of gas chromatography and mass spectrometry provides researchers with the ability to separate species before analyzing the mass spectrum.

Vernier Mini GC™ Plus and the IQ-2000™

The IQ-2000™ samples from the Mini GC™ Plus exhaust port, allowing students or instructors to quickly connect two analytical techniques, thereby providing in-depth compound identification.

The ability to easily connect or disconnect the instruments gives students the opportunity to fully distinguish the difference between each instrument’s capabilities, limitations, and sample requirements.

Using the retention time from the GC and the spectral data from the mass spectrometer, students can correctly identify individual chemicals or mixtures.

Sample Data for a mixture of Ethyl Acetate and Methyl Isobutyl Ketone

Top: Vernier Logger Pro® software displaying sample intensity of two organic solvents vs retention time.

Middle: Merlin Automation software depicting the mass spectrum of Ethyl Acetate (detail A).

Bottom: Merlin Automation software depicting the mass spectrum of Methyl Isobutyl Ketone (detail B).




Figure 1. Mixture Sample Data

Gas Chromatograph Vernier Mini GC™ Plus**Mass Spectrometer IQ-2000™**

Carrier Gas	Air
Column	Restek MXT-1 Stainless (11 m) Capillary Column
Column Temperature Ramp	Maximum of 10°C / minute Warm-up time: <10 minutes
Pressure Regulation	Software Controlled
Sample Introduction	Syringe Injection
Temperature Regulation	Software Controlled 30 - 160°C
Onboard Detector	MEMS Chemi – Capacitive with two levels of sensitivity
Software	Computer – Logger Pro®
Stand Alone	LabQuest® 2
Dimensions	10.8 cm x 13.3 cm x 19.1 cm

Compounds for Analysis	Any gas or vapor
Detection Range	100 ppb – 100%
Sample Introduction	Flow-by of GC exhaust to heated fused silica transfer line
Detector	Electron Multiplier
Software	Merlin Automation
Dimensions	Without Cart: 77 cm x 51 cm x 63 cm With Cart: 114 cm x 52 cm x 65 cm

**Vernier Mini GC™ Plus
Compounds for Analysis**

Mini GC Compound Type	Typical Compounds	Range of Acceptable Boiling Points
Alcohols	C ₁ -C ₆	50-175°C
Aldehydes	C ₁ -C ₈	50-175°C
Amides	C ₃ -C ₄	150-250°C
Polar-Substituted Aromatic Hydrocarbons	C ₆ -C ₁₂	100-150°C
Carboxylic Acids	C ₁ -C ₄	100-150°C
Esters	C ₂ -C ₁₀	30-200°C
Ethers	C ₂ -C ₈	30-110°C
Heterocycles	C ₅ -C ₈	100-150°C
Halogenated Hydrocarbons	C ₁ -C ₈	30-120°C
Ketones	C ₃ -C ₈	50-175°C
Nitriles	C ₂ -C ₅	50-120°C
Nitroalkane/aromatic	C ₁ -C ₆	100-210°C
Phosphonates	C ₃ -C ₉	180-210°C
Sulfides	C ₄ -C ₈	90-150°C



Figure 2. Extrel® IQ-2000™ Educational Mass Spectrometer with Vernier Mini GC™ Plus

